

There is Power Mobility. And then there's

POWERFUL MOBILITY



CityBug



There is Power Mobility. And then there's

POWERFUL MOBILITY



MAX SPEED
4.0
MPH



LOAD CAPACITY
250
LBS



BATTERY
12 Ah
SLA



RANGE
8
MILES

Model	S11E
Weight Capacity	250 Lbs (115kg)
Seat	Folding Mobility Seat, 16"
Tire Type	Solid Tires
Drive Wheel	22.8cm x 6.5cm(PU)
Front Caster (Wheel)	20cm x 5cm(PU)
Max Speed	4 Miles/Hr/ 6KM/H
Battery Specifications	12V 12Ah x 2pcs
Battery Range	Up to 8 Miles
Charger Type	2Amp, 100/240VAC Charger
Controller Type	S-Drive 45Amp

Motor Type	270W
Weight w/ Battery	67 Lbs / 30.8KG
Weight W/O Battery	50 Lbs / 22.4KG
Battery Weight	18 Lbs / 8.4KG
Turning Radius	37" (95cm)
Suspension	None
Length	32" (82.5cm)
Width	20" (52cm)
Height	34" - 38" (87cm~97cm)

Seat Width	16" (40cm)
Seat Height	21" (54cm)
Seat Depth	12" (33cm)
Back Height	9" (22cm)
Wheel Base	21" (53.5cm)
Ground Clearance	2" (5cm)
Climbing Angle	3 Degree

(*) SPECIFICATION PERFORMANCE IS DETERMINED IN THE LABORATORY SETTING BASED ON IDEAL CONDITIONS (MINIMAL FRICTION, NO WIND AND CONTROLLED TEMPERATURE AND PRESSURE WITH AN AVERAGE STANDARD WEIGHT). ACTUAL PERFORMANCE WILL BE DIFFERENT THAN RESULTS OBTAINED DURING LABORATORY TESTING.

Call us at 1-800-944-1997 | sales@evrider.com | www.evrider.com