

There is Power Mobility. And then there's

# POWERFUL MOBILITY



Transport M



There is Power Mobility. And then there's

# POWERFUL MOBILITY

## Transport M

Approved for airline travel, the S19M is equipped with an 11.5AH lithium battery and weighs just 42 lbs making it perfect for travel on planes, trains, and cruise ships. Due to its ultra-tight turning radius, this compact scooter offers great maneuverability in narrow spaces. It features a height-adjustable tiller, LED headlight, and dual hand two-way throttle control. It can be stored virtually anywhere and easily fits in the trunk of most automobiles.



MAX SPEED  
**3.8**  
MPH



LOAD CAPACITY  
**250**  
LBS



BATTERY  
**11.5 Ah**  
LITHIUM



RANGE  
**10**  
MILES

|                          |  |
|--------------------------|--|
| WEIGHT CAPACITY          | 250 lbs                                    |
| DRIVE WHEEL              | 8" x 2", Solid Tire, Flat Free             |
| FRONT WHEEL              | 7" x 1.6", Solid Tire, Flat Free           |
| REAR CASTER(ANTI-TIPPER) | yes, standard                              |
| MAX SPEED                | up to 3.8 mph *                            |
| BATTERY SPECIFICATIONS   | 25.2V 11.5Ah Lithium Battery x1 (included) |
| MAXIMUM TRAVEL RANGE     | up to 10 miles*                            |

|                 |  |
|-----------------|--|
| CHARGER TYPE    | 2 Amp 24 VDC, input: 110/240 Vac, Off Board Li-Ion Charger |
| CONTROLLER TYPE | Dynamic R-Series 50A                                       |
| MOTOR TYPE      | 180 watt 24 volt   |
| LENGTH          | 33.8"  |
| WIDTH           | 16.5"  |
| OVERALL HEIGHT  | 32.8"  |
| WEIGHT          | 42 lbs   |
| Colors          | Blue, Copper, Burgundy, Plum                               |

|                               |                       |
|-------------------------------|-----------------------|
| SUSPENSION                    | None                  |
| TURNING RADIUS                | 31"                   |
| SEAT HEIGHT                   | 20"                   |
| SEAT WIDTH                    | 16"                   |
| Seat Height – Cushion to Deck | 14"                   |
| WHEELBASE                     | 25"                   |
| GROUND CLEARANCE              | 1.4"                  |
| OVERALL DIMENSIONS            | 33.8" x 16.5" x 32.8" |
| FOLDED DIMENSIONS             | 28.7" x 17.3" x 16.5" |

(\*) SPECIFICATION PERFORMANCE IS DETERMINED IN THE LABORATORY SETTING BASED ON IDEAL CONDITIONS (MINIMAL FRICTION, NO WIND AND CONTROLLED TEMPERATURE AND PRESSURE WITH AN AVERAGE STANDARD WEIGHT). ACTUAL PERFORMANCE WILL BE DIFFERENT THAN RESULTS OBTAINED DURING LABORATORY TESTING.